

## Straight Laryngeal Mirror

A range of accessories make this instrument particularly versatile.



Metal tongue-depressor



Transilluminator



Protective cap for bulb



Laryngeal mirror	2.5 V	3.5 V
complete with plane mirror 20 mm dia .	<b>B-001.12.101</b>	<b>B-002.12.101</b>
same, without mirror	<b>B-001.12.109</b>	<b>B-002.12.109</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.057</b>	<b>X-002.88.063</b>

### Accessories

Plane mirror 20 mm dia.	<b>B-000.12.103</b>
Concave mirror 20 mm dia.	<b>B-000.12.104</b>
Metal tongue-depressor	<b>B-000.12.105</b>
Transilluminator	<b>B-000.12.106</b>
Protective cap for bulb (for mouth illumination)	<b>B-000.12.107</b>

## Curved Laryngeal Mirror

Assists examination with improved view of the larynx.



Metal tongue-depressor



Laryngeal mirror	2.5 V	3.5 V
complete with plane mirror 20 mm dia .	<b>B-001.12.201</b>	<b>B-002.12.201</b>
same, without mirror	<b>B-001.12.208</b>	<b>B-002.12.208</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>	<b>X-002.88.062</b>

### Accessories

Plane mirror 20 mm dia.	<b>B-000.12.203</b>
Concave mirror 20 mm dia.	<b>B-000.12.204</b>
Metal tongue-depressor	<b>B-000.12.205</b>

## HEINE mini 3000® Laryngeal Mirror



- Can only be used with the mini 3000 handle system.
- ∴ **A range of accessories make this instrument particularly versatile.**
- ∴ **Enhanced XHL Xenon Halogen Technology.**  
Concentrated illumination of oral cavity and pharynx.
- ∴ **Solid metal construction.** Durable and tough.

mini 3000 laryngeal mirror	2.5 V
with plane mirror 20 mm dia, without handle	<b>D-001.77.101</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.057</b>

## Tongue-Blade Holder

For wooden tongue-blades



- ∴ **XHL Xenon Halogen Technology.**  
Concentrated illumination of oral cavity and pharynx.
- ∴ **Light-guide for bulb.** Avoids dazzle.

<b>Tongue-Blade Holder</b>	2.5 V	3.5 V
for wooden tongue-blades	<b>B-001.12.302</b>	<b>B-002.12.302</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>	<b>X-002.88.049</b>
Replacement light-guide for bulb		<b>B-000.12.308</b>

## Tongue-Blade Holder

With Blade Ejector for use with HEINE Disposable Tongue-Blades



- ∴ **XHL Xenon Halogen Technology.**  
Concentrated illumination of oral cavity and pharynx.
- ∴ **Light-guide for bulb.** Avoids dazzle.
- ∴ **Blade ejector.** No need to touch used blades.
- ∴ **Only for HEINE Disposable Tongue-Blades.**

<b>Tongue-Blade Holder</b>	2.5 V	3.5 V
with blade ejector	<b>B-001.12.305</b>	<b>B-002.12.305</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>	<b>X-002.88.049</b>
Screw-on bulb-sleeve		<b>B-000.12.308</b>

### Disposable tongue-blades

Pack of 100 disposable tongue-blades	<b>B-000.12.304</b>
--------------------------------------	---------------------

- ∴ **Strong and hygienic.**





## HEINE alpha+® Tongue-Blade Holder

### With blade ejector for HEINE Disposable Tongue-Blades



- ∴ **XHL Xenon Halogen Technology.**  
Ideal illumination of mouth and throat.
- ∴ **Finger-tip pressure releases the used blade.**  
No need to touch the blade after use.
- ∴ **Only for HEINE Disposable Tongue-Blades.**

<b>alpha+ Tongue-Blade Holder</b>	2.5 V	3.5 V
with alpha+ battery handle	<b>D-001.74.120</b>	
same, with alpha+ rechargeable handle		<b>D-002.74.121</b>
same, with alpha+ rechargeable handle and alpha+ NT 200 B charger		<b>D-002.74.114</b>
alpha+ Tongue-Blade Holder without handle	<b>D-001.74.100</b>	<b>D-002.74.100</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>	<b>X-002.88.049</b>

#### Disposable tongue-blades

Pack of 100 disposable blades	<b>B-000.12.304</b>
-------------------------------	---------------------

**Strong and hygienic.**

## HEINE mini 3000® Tongue-Blade Holder

### With blade ejector for HEINE Disposable Tongue-Blades



- Can only be used with the mini 3000 handle system.
- ∴ **Enhanced XHL Xenon Halogen Technology.** Concentrated white light.  
Concentrated illumination of oral cavity and pharynx.
- ∴ **Attachment clip with integrated on/off switch.** Secure. Switches off automatically when replaced in the pocket.
- ∴ **Finger-tip pressure releases the used blade.**  
No need to touch the blade after use.
- ∴ **Solid metal construction.** Durable and tough.
- ∴ **Only for HEINE Disposable Tongue-Blades.**
- ∴ **2-piece handle and head.** Easy to maintain, flexible.
- ∴ **High-quality handle casing: Chrome-finish upper section/refined plastic.**  
Shockproof, sturdy, non-slip.
- ∴ **Replaceable batteries.** Size AA.

<b>mini 3000 Tongue-Blade Holder</b>	2.5 V
Complete with mini 3000 battery handle and 5 disposable blades	<b>D-001.74.118</b>
Tongue-blade holder head with bulb, without handle	<b>D-001.74.100</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.037</b>

#### Disposable tongue-blades

Pack of 100 disposable blades	<b>B-000.12.304</b>
-------------------------------	---------------------

**Strong and hygienic.**

HEINE mini 3000® pocket instruments – a complete range of quality diagnostic instruments for your pocket. Available in black or blue. Please specify the colour you prefer when ordering. If not otherwise specified, we will supply black instruments.

## Multi-purpose Illuminator



- **XHL Xenon Halogen Technology.** Bright, concentrated light.
- **Wide range of applications.** Used as a diagnostic light or tongue-depressor for wooden or disposable blades.

Multi-purpose Illuminator	2.5 V	3.5 V
Multi-purpose Illuminator	<b>B-001.12.320</b>	<b>B-002.12.320</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.077</b>	<b>X-002.88.078</b>

## Fiber Optic Nasal Illuminator



Accessory for Multi-purpose Illuminator.

- **Fiber Optic Nasal Illuminator.** For specula and spreadable nasal speculum.

Fiber Optic Nasal Illuminator	<b>B-000.12.321</b>
-------------------------------	---------------------

## Fiber Optic Nasal Illuminator, complete



Fiber Optic Nasal Illuminator	2.5 V	3.5 V
complete with nasal speculum	<b>B-001.12.322</b>	<b>B-002.12.322</b>
same, with spreadable nasal speculum	<b>B-001.12.323</b>	<b>B-002.12.323</b>
XHL Xenon Halogen spare bulb	<b>X-001.88.077</b>	<b>X-002.88.078</b>
Nasal speculum for above		<b>B-000.11.220</b>
Spreadable nasal speculum for above		<b>B-000.11.231</b>



HEINE ENT Diagnostic Sets



Set complete with **BETA 100 Diagnostic Otoscope**, 3 reusable specula, tongue-depressor for wooden tongue-blades and one spare bulb in hard case with:

	2.5 V	3.5 V
BETA battery handle	<b>B-242.10.118</b>	
BETA R rechargeable handle		<b>B-242.20.376</b>
BETA TR rechargeable handle and transformer		<b>B-242.20.384</b>
BETA NT rechargeable handle and NT 200 charger		<b>B-242.20.420</b>



Set complete with **BETA 100 Diagnostic Otoscope**, 1 set of reusable specula, 10 disposable specula, straight laryngeal mirror with accessories, spreadable nose speculum and one spare bulb each in zipper case with:

	2.5 V	3.5 V
BETA battery handle	<b>B-243.10.118</b>	
BETA R rechargeable handle		<b>B-243.20.376</b>
BETA TR rechargeable handle and transformer		<b>B-243.20.384</b>
BETA NT rechargeable handle and NT 200 charger		<b>B-243.20.420</b>



Set complete with **BETA 100 Diagnostic Otoscope**, 3 reusable specula, tongue-depressor for disposable tongue-blades, straight laryngeal mirror without accessories, spreadable nose speculum and one spare bulb each in hard case with:

	2.5 V	3.5 V
BETA battery handle	<b>B-244.10.118</b>	
BETA R rechargeable handle		<b>B-244.20.376</b>
BETA TR rechargeable handle and transformer		<b>B-244.20.384</b>
BETA NT rechargeable handle and NT 200 charger		<b>B-244.20.420</b>



Set complete with **K 100 Diagnostic Otoscope**, 3 reusable specula, tongue-depressor and one spare bulb for each instrument in hard case with:

	2.5 V	3.5 V
BETA battery handle	<b>B-241.10.118</b>	
BETA R rechargeable handle		<b>B-241.20.376</b>
BETA TR rechargeable handle with transformer		<b>B-241.20.384</b>
BETA NT rechargeable handle with NT 200 charger		<b>B-241.20.420</b>



Please note: BETA 3.5 V rechargeable handles can be supplied with a Lithium battery (**X-007.99.381**) with increased capacity and no memory effect at extra cost. Please specify when ordering.

To order sets with 3.5V BETA handle with Li-ion battery technology, use the number "27" in the voltage designation. E. g. B-242.**27**.420 for the BETA 100 Diagnostic Otoscope Set with BETA NT handle and charger.